# GEN<sup>X</sup>TRACT Blood DNA Extraction System

REF

2-014

IVD

100 extractions

2-8°C

The **GEN<sup>X</sup>TRACT Blood** DNA Extraction System is designed to be used for convenient and rapid extraction of DNA from anticoagulated (EDTA/citrate) blood as template for *in vitro* amplification (PCR)

- 
- |                                      |          |
|--------------------------------------|----------|
| 1. <b>Lysis Solution</b>             | 4x 50 ml |
| 2. <b>GEN<sup>X</sup>TRACT Resin</b> | 4x 5 ml  |

*Resuspend each time immediately before removing an aliquot.*

---

## MATERIALS REQUIRED BUT NOT SUPPLIED

- Microtubes (1.5 ml with screw cap)
- Adjustable micropipettes
- Adjustable microcentrifuge capable of 3,000-12,000 rpm (approx. 1,000-12,000 x g)
- Incubator (e.g. heating block, water bath) capable of 56°C and 98°C (± 2°C)

### ViennaLab Diagnostics GmbH

Gaudenzdorfer Guertel 43-45

A-1120 Vienna, Austria

Phone: (+43-1) 8120156-0

Fax: (+43-1) 8120156-19

info@viennalab.co.at

ESTABLISHED INNOVATIONS IN DIAGNOSTICS

www.viennalab.com

## ASSAY PROCEDURE

Use fresh or frozen blood with EDTA or citrate anticoagulant; avoid blood containing heparin. Do not store blood for more than 3 days at ambient temperature or more than 1 week at 2-8°C before use. Blood which has been kept frozen for more than one year, or gone through more than three freeze-thaw cycles is unsuitable to be used in this procedure. Bring blood samples to room temperature. Mix well by carefully inverting blood collection tubes several times. Repeat mixing each time before withdrawing an aliquot of blood. Allow Lysis Solution and GEN<sup>X</sup>TRACT Resin to reach room temperature.

- Pipette **100 µl blood sample** into a 1.5 ml microtube with screw cap.
- Add **1 ml Lysis Solution**, close tube and mix by inverting several times.
- Let stand for **15 min.** at room temperature.
- Centrifuge for **5 min.** at **3,000 rpm** (approx. 1,000 x g) in a microcentrifuge.
- Remove and discard the upper (top) 1 ml of supernatant.
- Add **1 ml Lysis Solution**, close tube and mix by inverting several times.
- Centrifuge for **5 min.** at **12,000 rpm** (approx. 12,000 x g) in a microcentrifuge.
- Remove and discard the supernatant except for approx. 50 µl of a visible, soft pellet.
- Resuspend GEN<sup>X</sup>TRACT Resin by swirling the bottle thoroughly.
- Add **200 µl GEN<sup>X</sup>TRACT Resin** to the pellet. Close tube and vortex for 10 sec.  
⚠ GEN<sup>X</sup>TRACT Resin sediments quickly. Repeat resuspension each time immediately before removing another aliquot.
- Incubate for **20 min.** at **56°C**. Vortex for 10 sec.
- Incubate for **10 min.** at **98°C**. Vortex for 10 sec.
- Centrifuge for **5 min.** at **12,000 rpm** in a microcentrifuge. Cool on ice.

The resulting supernatant contains DNA template suitable for immediate use in PCR. For further storage, the supernatant should be transferred into a fresh tube and kept refrigerated (2-8°C; up to one week) or frozen at -20°C.

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