

# NomeIRT™

## Western Blot Stripping Buffer

RUO Research Use Only

REF 21112



### INTRODUCTION

Western blotting is widely used to detect and compare proteins in complex mixtures, and chemiluminescence has largely replaced chromogenic substrates as the most convenient and sensitive method of detection. Nitrocellulose and PVDF membranes probed by Western blotting procedures and detected by chemiluminescent or other non-precipitating substrates can be stripped and re-probed using NomeIRT™ Western Blot Stripping Buffer. One advantage of chemiluminescence is the ability to strip and reprobe the protein mixture on the membrane. Traditional stripping methods use conditions that are effective for only low-affinity antibody-antigen interactions or are so harsh that they tend to adversely alter the antigen for subsequent immunoprobings. No stench smelling and without addition of 2-mercaptoethanol or its analogs. Incubate the blot in NomeIRT™ Western Blot stripping reagent at room temperature.

### KIT CONTENTS

Label	Contain
NomeIRT™ Western Blot Stripping Buffer	500 ml

### STORAGE AND STABILITY

- Storage condition : Upon arrival store this product at 4°C ~ room temperature. For long term, store the product at 4°C. Product shipped at ambient temperature.

### APPLICATIONS

- Western Blot

### ADDITIONAL REQUIRED EQUIPMENT

- Shakers
- Shaking incubator (optional)

### PROTOCOL

After initial probing, be sure to keep membrane wet in TBST buffer. DO NOT DRY!

- Warm the bottle of NomeIRT™ Western Blot Stripping Buffer to room temperature.
- Pour 15~30 ml stripping reagent to a clean container and put the blot in the container. Make sure that the blot is fully submerged with the stripping buffer.
- Incubate the blot in stripping reagent at room temperature for 10 ~ 30 minutes with agitation. Though incubation with the high affinity antibodies need to be optimized, 15 minutes stripping at room temperature is usually sufficient for most of antibodies.

**Note:** Optimization of both incubation time and temperature is essential for best results. In general, higher affinity antibodies will require at least 30 minutes of stripping and may require an incubation temperature of 37°C.

- Wash for 5 minutes in TBS-T at room temperature using large volumes (e.g. 100 ml) of wash buffer.

**Note:** To test the stripping effect, pour ECL reagent on blot followed by 5 minutes exposure to a film.

### ORDERING INFORMATION

Product Name	Amount	Cat. No.
GangNam-STAIN™ Prestained Protein Ladder	250 ul	24052
WEST-lott™ PVDF Transfer membrane (0.45 um, 0.22 um)	300 x 3000 mm	ITM-P3031 / ITM-P3032
Ponceau Ssolution	1L	31192
WEST-ZOL plus Western Blot Detection System	200 ml	16024

#### ◆ References

- Kaufmann, S.H., et al. (1987). The erasable Western blot. *Anal Biochem* 161: 89-95.
- Kaufmann, S.H. and Kellner, U. (1998). Erasure of Western blots after autoradiographic or chemiluminescent detection. In *Immunocytochemical Protocols*. Ed. POUND, J.D. Humana Press, Totowa, N.J., p. 223-35.

Technical support : +82-505-550-5600

Review date : 2014.06.26