α -Globin StripAssay[®]

REF

4-160

IVD

10 Tests

2-8°C

CE

1a.	Amplification Mix A1 (<i>yellow cap</i>)	250 μ l	
1b.	Amplification Mix A2 (<i>white cap</i>)	250 μ l	
1c.	Amplification Mix B (<i>green cap</i>)	250 μ l	
2.	Taq Dilution Buffer (<i>transparent cap</i>)	500 μ l	
3.	HS-Taq DNA Polymerase (5 U/μl) (<i>red cap</i>)	125 U	
4.	DNAT (<i>blue cap</i>)	1.5 ml	Warning
5.	Typing Trays	3	
6a.	Teststrips A (<i>black cap</i>)	10	
6b.	Teststrips B (<i>white cap</i>)	10	
7.	Hybridization Buffer (<i>white cap</i>)	25 ml	
8.	Wash Solution A (<i>white cap</i>)	80 ml	
9.	Conjugate Solution	25 ml	
10.	Wash Solution B	80 ml	
11.	Color Developer	25 ml	Warning
12.	Instructions For Use	1	
13.	Collector[™] Sheet	1	

ViennaLab Diagnostics GmbH

Gaudenzdorfer Guertel 43-45

A-1120 Vienna, Austria

Phone: (+43-1) 8120156-0

info@viennalab.com

www.viennalab.com


Fig. 1: Teststrip Design

Note: Teststrips are not drawn in real size and must not be used for interpretation of results!

Instructions for use

I. INTENDED USE

Assay for the identification of α -globin gene mutations based on polymerase chain reaction (PCR) and reverse-hybridization. *For human in vitro diagnostics.*

II. METHODOLOGY

The procedure includes three steps: (1) DNA isolation, (2) PCR amplification using biotinylated primers, (3) hybridization of amplification products to test strips containing allele-specific oligonucleotide probes immobilized as an array of parallel lines (Fig. 1). Bound biotinylated sequences are detected using streptavidin-alkaline phosphatase and color substrates.

The assay covers 21 α -globin mutations: 3.7 single gene deletion, 4.2 single gene deletion, MED-I double gene deletion, SEA double gene deletion, THAI double gene deletion, FIL double gene deletion, 20.5 kb double gene deletion, anti-3.7 gene triplication, $\alpha 1$ cd 14 [TGG>TAG], $\alpha 1$ cd 59 [GGC>GAC] (Hb Adana), $\alpha 2$ init cd [ATG>ACG], $\alpha 2$ cd 19 [-G], $\alpha 2$ IVS1 [-5nt], $\alpha 2$ cd 59 [GGC>GAC], $\alpha 2$ cd 125 [CTG>CCG] (Hb Quong Sze), $\alpha 2$ cd 142 [TAA>CAA] (Hb Constant Spring), $\alpha 2$ cd 142 [TAA>AAA] (Hb Icaria), $\alpha 2$ cd 142 [TAA>TAT] (Hb Pakse), $\alpha 2$ cd 142 [TAA>TCA] (Hb Koya Dora), $\alpha 2$ poly A-1 [AATAAA-AATAAG], $\alpha 2$ poly A-2 [AATAAA-AATGAA]

Further genetic information is available at OMIM Online Mendelian Inheritance in Man: www.ncbi.nlm.nih.gov/omim

III. KIT COMPONENTS

See list of all kit components on page I.

Warning: DNAT contains 1.6% NaOH.

H315: Causes skin irritation

H319: Causes serious eye irritation

P280: Wear protective gloves/protective clothing/eye protection/face protection

P337 + P313: If eye irritation persists: Get medical advice/attention

Warning: Color Developer contains $\leq 0.4\%$ maleic acid.

H317: May cause an allergic skin reaction

P280: Wear protective gloves/protective clothing/eye protection/face protection

P302 + P352: If on skin: wash with plenty of water

P333 + P313: If skin irritation or rash occurs: get medical advice /attention

Amplification Mix, Taq Dilution Buffer, Conjugate Solution, Wash Solution B contain 0.05% NaN_3 . Conjugate Solution contains streptavidin-alkaline phosphatase. Color Developer contains nitro blue tetrazolium (NBT) and 5-bromo-4-chloro-3-indolyl phosphate (BCIP).

Store all reagents at 2-8°C when not in use !

IV. MATERIALS REQUIRED BUT NOT SUPPLIED

In addition to standard molecular biology laboratory equipment, the following is needed:

- Thermocycler and suitable thin-walled plastic reaction tubes/strips
- Waterbath with shaking platform and adjustable temperature ($45^\circ\text{C} \pm 1^\circ\text{C}$)
- Vacuum aspiration apparatus
- Shaker (rocker or orbital shaker)
- *Optional: agarose gel electrophoresis equipment (for control of amplification products)*

V. ASSAY PROCEDURE**1. DNA Isolation**

It is recommended to use the following kit for DNA isolation from whole blood and various other sources:

REF 2-020: Spin Micro DNA Extraction Kit

The use of other DNA extraction methods with the α-Globin StripAssay® has not been evaluated.

2. In Vitro Amplification (PCR; 3 separate reactions per sample)

Keep all PCR reagents and DNA templates refrigerated throughout. Perform all steps until start of the thermal cycling program on ice (0-4°C).

- Prepare a fresh working dilution (1:15, final conc. 0.33 U/μl) of **HS-Taq DNA Polymerase** (red cap) in **Taq Dilution Buffer** (transparent cap).
- Prepare three reaction tubes for each sample to be amplified. Place tubes on ice.
- For each sample prepare 3 final PCR reaction mixes (A1, A2, B) on ice:
 - A1: **15 μl Amplification Mix A1** (yellow cap)
5 μl diluted HS-Taq DNA Polymerase (1.66 U)
5 μl DNA template
 - A2: **15 μl Amplification Mix A2** (white cap)
5 μl diluted HS-Taq DNA Polymerase (1.66 U)
5 μl DNA template
 - B: **15 μl Amplification Mix B** (green cap)
5 μl diluted HS-Taq DNA Polymerase (1.66 U)
5 μl DNA template

If DNA templates not prepared by the Spin Micro DNA Extraction Kit are used, a DNA concentration range of 2-10 ng/μl (= 10-50 ng DNA per reaction) is recommended.

- *Cap tubes tightly. Preheat the thermocycler to 95°C.*
- Insert reaction tubes and run the following thermocycling program:
 - pre-PCR: 95°C/5 min.**
 - thermocycling: 97°C/40 sec. - 64°C/40 sec. - 72°C/1:30 min. (3 cycles)**
97°C/40 sec. - 58°C/40 sec. - 72°C/1:30 min. (37 cycles)
 - final extension: 72°C/5 min.**

Store amplification products on ice or at 2-8°C for further use.

Optional: *Analyze amplification products by gel electrophoresis (e.g. 3% agarose gel).*

Fragment lengths: 881 bp; deletions: 1783 bp (A1)

296 bp; deletions: 250, 329, 373, 474, 578, 1025 bp (A2)

302, 864 bp; deletions: 1772 bp (B)

3. Hybridization - 2 Teststrips per sample (45°C; shaking waterbath)

Adjust the water level of the waterbath to approx. ½ of the height of the Typing Tray.

Heat the waterbath to exactly 45°C ($\pm 1^\circ\text{C}$). Check water temperature with a calibrated thermometer.

Prewarm Hybridization Buffer and Wash Solution A to 45°C. (Take care that all precipitates formed at 2-8°C become completely dissolved.)

Allow Teststrips, DNAT, Conjugate Solution, Wash Solution B and Color Developer to reach room temperature. Prepare Typing Tray(s).

Remove one Teststrip A and one Teststrip B for each sample using clean tweezers. (Touch Teststrips with gloves only!) Label Teststrips outside of the marker lines with a pencil. (No ballpoint pens, markers, etc.)

For all **Teststrips A** (one lane per sample):

- Pipette **20 µl DNAT** (blue cap) into the lower corner of each lane to be used in the Typing Trays.
- Add **10 µl amplification product A1** into the corresponding drop of DNAT.
Add **10 µl amplification product A2** into the same drop.
Mix thoroughly with a pipette. (The solution will remain blue.)

For all **Teststrips B** (one lane per sample):

- Pipette **10 µl DNAT** (blue cap) into the lower corner of each lane to be used in the Typing Trays.
- Add **10 µl amplification product B** into the corresponding drop of DNAT.
Mix thoroughly with a pipette. (The solution will remain blue.)

- Let stand for **5 min.** at room temperature.
- Add **1 ml Hybridization Buffer** (prewarmed to 45°C) into each lane.
Gently agitate tray. (The blue color will disappear.)
- Insert **Teststrip A** or **Teststrip B** with marked side up (lines visible!) into the respective lanes. Submerge completely.
- Incubate for **30 min.** at **45°C** on the shaking platform of the waterbath.
Set moderate shaking frequency (approx. 50 rpm) to avoid spilling. Keep the cover of the waterbath closed to avoid variations in temperature.
- At the end of incubation remove hybridization solutions by vacuum aspiration.
Proceed immediately. Do not allow Teststrips to run dry during the entire procedure.

4. Stringent Wash (45°C; shaking waterbath)

- Add **1 ml Wash Solution A** (prewarmed to 45°C). Rinse briefly (10 sec.).
Remove liquids by vacuum aspiration.
- Add **1 ml Wash Solution A** (45°C).
- Incubate for **15 min.** at **45°C** in the shaking waterbath.
Remove liquids by vacuum aspiration.
- Add **1 ml Wash Solution A** (45°C).
- Incubate for **15 min.** at **45°C** in the shaking waterbath.
Remove liquids by vacuum aspiration.

5. Color development (room temperature)

- Add **1 ml Conjugate Solution**.
- Incubate for **15 min.** at **room temperature** on a rocker or orbital shaker. Remove liquids by vacuum aspiration.
- Add **1 ml Wash Solution B**. Rinse briefly (10 sec.). Remove liquids by vacuum aspiration.
- Add **1 ml Wash Solution B**.
- Incubate for **5 min.** at **room temperature** on a rocker or orbital shaker. Remove liquids by vacuum aspiration.
- Add **1 ml Wash Solution B**.
- Incubate for **5 min.** at **room temperature** on a rocker or orbital shaker. Remove liquids by vacuum aspiration.
- Add **1 ml Color Developer**.
- Incubate for **15 min.** at **room temperature in the dark** on a rocker or orbital shaker. *A purple staining will appear upon positive reaction.*
- Wash Teststrips several times with distilled water. Let strips dry in the dark on absorbent paper. *Do not expose Teststrips to intense light after Color Development.*

VI. INTERPRETATION OF RESULTS

The genotype of a sample is determined from corresponding Teststrips A and B using the enclosed Collector™ sheet.


Place both processed Teststrips into the designated fields, align them to the schematic drawing using the red marker line (top) and the green or blue marker line (bottom), and fix them with adhesive tape.

A positive reaction of the uppermost Control line indicates the correct function of Conjugate Solution and Color Developer. This line should always stain positive.

A positive reaction of the PCR Control A and PCR Control B lines indicates the presence of the correct amplification products. These lines should always stain positive, except for the "negative control" which contains water in place of DNA template (see example H, page III).

For each polymorphic position, one of the following staining patterns should be obtained:

Note: Staining intensities of positive lines may vary. This is of no significance for the result.

	wild type line	mutant line	genotype
NOR	positive	negative	normal
HET	positive	positive	heterozygous
HOM	negative	positive	homozygous mutant

See examples of StripAssay results on page III (Fig. 3).

Some of the point mutations covered by the α-Globin StripAssay® are located within a few nucleotides on the α-globin gene. On the Teststrips these are represented by a common wild type probe, so that the 21 mutations are covered by 9 wild type probes only.

line	wild type probe	mutation
10	α1 cd 14	α1 cd 14 [G>A]
11	α1 cd 59	Hb Adana
25	α2 init cd	α2 init cd [T>C]
26	α2 cd 19	α2 cd 19 [-G]
27	α2 IVS 1	α2 IVS 1 -5nt
28	α2 cd 59	α2 cd 59 [G>A]
29	α2 cd 125	Hb Quong Sze
30	α2 cd 142	Hb Constant Spring, Hb Icaria, Hb Pakse, Hb Koya Dora
31	α2 poly A	α2 poly A-1, α2 poly A-2

Samples that are compound heterozygous for two of these mutations (e.g. Hb Constant Spring + Hb Pakse) will be lacking the common wild type signal (see example E, page III).

Samples that are compound heterozygous for one of the α1/α2 mutations, and a single or double gene deletion, will in most cases be lacking the respective wild type signal (see example D, page III).

For single and double gene deletions, several wild type probes discriminate between the heterozygous and the homozygous mutant state (see examples B and C, page III).

deletion	heterozygous	homozygous mutant
- 3.7	all WT signals present	WT signals 25-31 absent
- 4.2	all WT signals present	WT signals 25-31 absent
- 20.5 kb	all WT signals present	WT signals 10 and 25-31 absent
-- MED-I	all WT signals present	all WT signals absent
-- SEA	all WT signals present	all WT signals absent
-- THAI	all WT signals present	all WT signals absent
-- FIL	all WT signals present	all WT signals absent

The α-Globin StripAssay® does not allow to distinguish between the heterozygous and the homozygous mutant state of the anti-3.7 gene triplication (anti-3.7/αα and anti-3.7/anti-3.7).

Advice on troubleshooting may be obtained by contacting ViennaLab through the local distributor or directly at techhelp@viennalab.com.

VII. QUALITY CONSIDERATIONS

- A thorough understanding of the procedure outlined here, and precise laboratory equipment and techniques are required to obtain reliable results. Use of the StripAssay for human *in vitro* diagnostics needs to be limited to appropriately trained personnel.
- Do not use StripAssay components beyond the expiration date printed on the outside of the kit box. Do not mix reagents from different lots.
- Avoid microbial contamination and cross-contamination of reagents or samples by using sterile disposable pipette tips throughout. Do not interchange bottle caps.
- The Control line immobilized on each Teststrip allows a performance control of the chromogenic detection system. To monitor and validate the specificity of the hybridization and washing steps, control DNAs of known genotype should be included into each individual experiment.

VIII. SAFETY

- Do not drink, eat, smoke, or apply cosmetics in designated work areas. Wear laboratory coats and disposable gloves when handling specimens and kit reagents. Wash hands thoroughly afterwards.
- Handle specimens as if capable of transmitting infectious agents. Thoroughly clean and disinfect all materials and surfaces that have been in contact with specimens. Discard all waste associated with clinical specimens in a biohazard waste container.
- Avoid contact of DNAT with skin, eyes, or mucous membranes. If contact does occur, immediately wash with large amounts of water. If spilled, dilute with water before wiping dry.
- Adhere to all local and federal safety and environmental regulations which may apply.

IX. ASSAY PERFORMANCE SPECIFICATIONS

Performance data on 330 tested samples (660 alleles: 300 WT, 360 MUT).

Reference methods: DNA sequencing, gap PCR, multiplex gap-PCR, RDB, ARMS-PCR.

Sensitivity: 99.7% (95% CI: 98.46% to 99.95%)

(= probability of a positive test result in the presence of the examined sequence variant)

Specificity: 100% (95% CI: 98.77% to 100%)

(= probability of a negative test result in the absence of the examined sequence variant)

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Le test comprend 21 mutations de α-globine: 3.7 single gene deletion, 4.2 single gene deletion, MED-I double gene deletion, SEA double gene deletion, THAI double gene deletion, FIL double gene deletion, 20.5 kb double gene deletion, anti-3.7 gene triplication, α1 cd 14 [TGG>TAG], α1 cd 59 [GGC>GAC] (Hb Adana), α2 init cd [ATG>ACG], α2 cd 19 [-G], α2 IVS1 [-5nt], α2 cd 59 [GGC>GAC], α2 cd 125 [CTG>CCG] (Hb Quong Sze), α2 cd 142 [TAA>CAA] (Hb Constant Spring), α2 cd 142 [TAA>AAA] (Hb Icaria), α2 cd 142 [TAA>TAT] (Hb Pakse), α2 cd 142 [TAA>TCA] (Hb Koya Dora), α2 poly A-1 [AATAAA-AATAAG], α2 poly A-2 [AATAAA-AATGAA].

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