# IL28B RealFast™ Assay

REF 7-200 / 7-203 100 / 32 reactions  
-30°C -15°C CE IVD

**ViennaLab Diagnostics GmbH**  
Gaudenzdorfer Guertel 43-45  
A-1120 Vienna, Austria  
Phone: (+43-1) 8120156-0  
[info@viennalab.com](mailto:info@viennalab.com)  
[www.viennalab.com](http://www.viennalab.com)

## 1. Intended Use

The IL28B RealFast™ Assay is a fast and accurate real-time PCR test for the detection of the TT>ΔG dinucleotide frame-shift variant in the human *interferon lambda 4 (IFNL4)* gene/pseudogene, which is located upstream of the gene encoding interleukin 28B (IL28B). Genotyping helps to predict the therapeutic response in Hepatitis C Virus (HCV) infected patients who are being considered for pegylated-interferon/ribavirin (pegIFN/RBV) therapy. Homozygous carriers of the favorable TT allele are predisposed to an up to threefold higher sustained virological response (SVR) than carriers of the ΔG allele. In a human DNA extract the qualitative assay discriminates the three possible IFNL4 TT>ΔG genotypes associated with high (TT/TT) or low (TT/ΔG and ΔG/ΔG) treatment-induced or spontaneous HCV clearance. Reference sequence: HGVS: NG\_042193.1:g.1457\_1458delAAinsC; NCBI dbSNP: rs368234815.

## 2. Introduction

Chronic HCV infection affects around 3% of the world population and is the leading cause of cirrhosis, liver cancer and liver transplantation. Success of the recommended standard pegIFN/RBV therapy depends on both viral and host factors. Patients infected with the unfavorable HCV type 1 or 4 have lower SVR than those infected with type 2 or 3. Host factors influencing treatment response include age, gender, ethnicity and genetic polymorphisms affecting antiviral activity. The *IFNL4* rs368234815 TT>ΔG variant is strongly linked to rs12979860 C>T, a well-established marker for long-term therapeutic success. Several studies demonstrated rs368234815 as the functional variant for IL28B mRNA regulation and, compared to rs12979860, as the better predictor of SVR and HCV clearance particularly in patients with African ancestry. The IFNL4 genotype TT/TT favors spontaneous as well as treatment-induced viral clearance in HCV patients, primarily in those infected with type 1 or 4.

## 3. Kit Contents

100 / 32 Rxn

RealFast™ 2x Genotyping Mix	1 vial	□ white cap	1000 / 320 µl
IL28B Assay Mix	1 vial	■ purple cap	550 / 550 µl
IL28B TT/TT-Control	1 vial	■ green cap	75 / 75 µl
IL28B ΔG/ΔG-Control	1 vial	■ red cap	75 / 75 µl

The RealFast™ 2x Genotyping Mix comprises HotStart Taq DNA polymerase and dNTPs in an optimized buffer system.

The IL28B Assay Mix consists of *IFNL4* gene-specific primers and two allele-specific, dual-labeled hydrolysis probes. Controls representing IL28B TT/TT and ΔG/ΔG genotypes are supplied with the kit.

The kit contains reagents for 100 / 32 reactions in a final volume of 20 µl each.

## 4. Storage and Stability

IL28B RealFast™ Assay is shipped on cooling blocks. On arrival, store the kit at -30 to -15°C. Alternatively, store at 2 to 8°C for short-term use within one month. The kit withstands up to 20 freeze/thaw cycles with no loss of activity. Avoid prolonged exposure to intense light. If stored correctly, the kit will retain full activity until the expiration date indicated on the label.

## 5. Product Description

### 5.1. Principle of the Test

The test is based on the fluorogenic 5' nuclease assay, also known as TaqMan® assay. Each reaction contains a gene-specific primer pair which amplifies a 110 bp fragment of the *IFNL4* gene, and two dual-labeled, allele-specific hydrolysis probes which hybridize to the target sequence of the amplified fragment. The proximity of the 5'-fluorescent reporter and 3'-quencher dye on intact probes prevents the reporter from fluorescing. During the extension phase of PCR the 5' – 3' exonuclease activity of the Taq DNA polymerase cleaves the 5'-fluorescent reporter from the hybridized probe. The physical separation of the fluorophore from the quencher dye generates a fluorescent signal in real-time, which is proportional to the accumulated PCR product.

In homozygous TT/TT samples the **HEX-labeled TT probe** hybridizes to the complementary strand of the gene fragment. A strong fluorescence signal is detected in the HEX channel (556nm) and no or only a baseline signal in the FAM channel (520nm). Vice versa, in homozygous ΔG/ΔG samples the **FAM-labeled ΔG probe** binds to the gene fragment. A strong fluorescence signal is detected in the FAM channel and no or only a baseline signal in the HEX channel. In heterozygous samples (TT/ΔG) both probes bind to the amplicons and generate intermediate signals in both channels.

### 5.2. Real-time PCR Instrument Compatibility

The IL28B RealFast™ Assay is validated for use with the AB 7500 Fast instrument.

The kit is compatible with various common real-time PCR instruments capable of recording FAM and HEX fluorescence:

- ✓ AB 7500 Fast (Applied Biosystems®)
- ✓ AB StepOne™ (Applied Biosystems®)
- ✓ CFX96™ (Bio-Rad)
- ✓ LightCycler® 480 (Roche)
- ✓ MIC qPCR Cyclers (bms)
- ✓ Mx3005P (Agilent Technologies)
- ✓ Rotor-Gene® 6000 (Qiagen)

» **Note:** RealFast™ Genotyping QuickGuides for setting up and analyzing experiments on different types of instruments can be downloaded from [www.viennalab.com](http://www.viennalab.com).

When using AB StepOne™, set passive reference dye to "ROX"! «

The kit is supplied **without ROX**. For use with real-time PCR instruments requiring high ROX for normalization of data (e.g. Applied Biosystems® instruments: StepOne™, 7300, 7900/7900HT), add ROX at a final concentration of 1 µM to the 2x Genotyping Mix.

### 5.3. Assay Performance Specifications

The IL28B RealFast™ Assay and a reference kit for the detection of the highly linked rs12979860 C>T polymorphism showed 99,4% concordance of genotyping results obtained on 84 samples.

Out of 168 alleles, the IL28B RealFast™ Assay identified 113 favorable TT and 55 unfavorable ΔG alleles, whereas the reference kit identified 114 favorable C alleles and 54 unfavorable T alleles. Sequencing of one discrepant sample revealed that the ΔG allele was present, but not correlated with the rs12979860 T allele on one chromosome.

The **sensitivity** and **specificity** of the IL28B RealFast™ Assay showed a 100% true positive and a 100% true negative rate, respectively.

Limit of detection: 0.2 ng genomic DNA (per reaction). Recommended DNA concentration: 2 to 20 ng/µl genomic DNA.

## 6. Materials Required but not Supplied

Real-time PCR instrument with FAM (520 nm) and HEX (556 nm) filters, instrument-compatible reaction vessels, disposable powder-free gloves, vortexer, mini-centrifuge for 2.0 ml tubes, tube racks, set of calibrated micropipettes (0.5 – 1000 µl), sterile tips with aerosol-barrier filter, molecular grade water, DNA extraction system, freezer, biohazard waste container.

## 7. Experimental Protocol

### 7.1. DNA Extraction

DNA extraction reagents are **not supplied** with the kit.

DNA isolated from various specimens (e.g. whole peripheral blood, dried blood spots, buccal swabs or saliva) can be used. Ensure extracted DNA is suitable for amplification in terms of concentration, purity and integrity.

For accurate genotype calling, the DNA amount per reaction should be within the range of 10 to 100 ng for all samples.

### 7.2. PCR Controls

**Always** include a **No Template Control (NTC)** in each experiment to confirm absence of potential contamination. It is advisable to run the NTC (use PCR-grade water instead of DNA) in duplicate.

**Always** include the IL28B **TT/TT-Control** and IL28B **ΔG/ΔG-Control** as positive reference signals for your unknown samples. Some real-time PCR software, e.g. AB 7500 Fast, requires signals for all three possible genotypes for correct allelic discrimination. In order to obtain a heterozygous control (TT/ΔG-Control), mix an aliquot of TT/TT-Control and ΔG/ΔG-Control in a ratio of 1:1.

» **Note:** TT/TT- and ΔG/ΔG-Controls are potential sources of contamination. Make sure to handle them carefully. «

### 7.3. Preparation of IL28B RealFast™ Master Mix

Gently vortex and briefly centrifuge all solutions after thawing. Set up PCR at room temperature. Prepare sufficient **Master Mix** for all your reactions (N samples + positive controls + negative controls) plus at least one additional reaction to compensate for pipetting inaccuracies:

Component	per reaction	e.g. 24+1 reactions
RealFast™ 2x Genotyping Mix	10 µl	250 µl
IL28B Assay Mix	5 µl	125 µl
<b>Master Mix</b>	<b>15 µl</b>	<b>375 µl</b>

Dispense **15 µl Master Mix** into each well. Add **5 µl purified DNA** or **Control** template to reach a final reaction volume of 20 µl.

To minimize risk of contamination, always pipette templates in the following order: first NTC, then samples, last positive controls. Immediately close reaction vessels.

» **Note:** Avoid creating bubbles in the final reaction mix and avoid touching the optical surface of the cap or sealing film without gloves. Both may interfere with fluorescence measurements. Centrifuge briefly if needed. «

### 7.4. PCR Program

Program the real-time PCR instrument according to the manufacturer's instructions for allelic discrimination / genotyping experiments. Place the samples into the thermal cycler and run the following program:

**AB 7500 Fast, StepOne™, CFX96™, LightCycler® 480, Mx3005P** and other Peltier heating block-based instruments:

Cycles	Temp	Time	Steps
1	95 °C	3 min	Initial denaturation
40	95 °C	15 sec	Denaturation
	<b>60 °C</b>	1 min	Annealing/Extension – <b>Data acquisition</b> on FAM and HEX channel

**MIC qPCR Cycler, Rotor-Gene® 6000\*):**

Cycles	Temp	Time	Steps
1	95°C	3 min	Initial denaturation
40	95°C	15 sec	Denaturation
	<b>60°C</b> *for 36-well rotor: <b>56°C</b>	1 min	Annealing/Extension – <b>Data acquisition</b> on Green and Yellow channel

## 8. Data Analysis / Interpretation of Results

The genotype of each sample is determined by calculating the ratio between signals recorded in the **HEX channel (TT)** and signals recorded in the **FAM channel (ΔG)**. Most real-time PCR software automatically resolves data of both channels into clusters in a scatterplot. Data points plotted along the x- and y-axes correspond to homozygous TT/TT and ΔG/ΔG genotypes, respectively. Data points clustered in the middle of the scatterplot represent heterozygous TT/ΔG genotypes. The NTC appears in the lower left corner.

Controls	Amplification in <b>FAM</b> channel (520 nm)	Amplification in <b>HEX</b> channel (556 nm)	Genotype / HCV clearance
TT/TT-Control	NO	<b>YES</b>	homozygous TT / high
TT/ΔG -Control	<b>YES</b>	<b>YES</b>	heterozygous TT/ΔG / low
ΔG/ΔG-Control	<b>YES</b>	NO	homozygous ΔG / low
NTC	NO	NO	----


Some instrument software needs manual threshold settings for accurate genotype calling.

Recommendations for Threshold Settings (C<sub>q</sub>):

Set threshold value for the FAM channel just above the background fluorescent signal generated by the TT/TT-Control (HEX-positive). Vice versa, set threshold value for the HEX channel just above the background fluorescent signal of the ΔG/ΔG-Control (FAM-positive).

Samples crossing the threshold line beyond C<sub>q</sub> 37 give invalid results and must be repeated.

To analyze acquired data, please follow your instrument software instructions.

## 9. Warnings and Precautions

- For *in vitro* diagnostics use only.
- Always use disposable powder-free gloves and wear suitable lab coat when handling specimens and reagents.
- Perform reaction setup in an area separate from nucleic acid preparation and PCR product analysis.
- Use pipettes dedicated for PCR setup only, use aerosol-guarded pipette tips.
- Use instrument-compatible reaction vessels with optically clear caps or sealers.
- Do not mix reagents from different lots.
- Do not use expired kits or kit components.

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