

Whitepaper

Rare Pathogens – Culture-Independent Eubacterial 16S rRNA Gene Diagnosis

Michael G. Lorenz – Molzym GmbH & Co. KG, Bremen, Germany

Keywords: Sepsis, septicaemia, bacteraemia, diversity, Gram-positive, Gram-negative pathogens, direct diagnosis, ribosomal genes, broad-range PCR, unbiased Real-Time PCR, Sanger sequencing, routine diagnosis, fluid and tissue biopsies, one sample processing protocol, contamination-free products, automated extraction

Summary

Rare pathogens are associated with critical conditions of patients in many cases demanding timely identification for adjustment of the antibiotic regime. Cultural diagnosis of rare pathogens, however, can be laborious and time-consuming, which lowers the prognosis of successful treatment. A prominent reason for this is the enormous diversity of rare pathogens that may be identifiable only with delay or even not identifiable by usual phenotyping methods. A review of clinical studies conducted in five countries from Europe, Asia and the USA reveals that 90.7 to 98.2% of systemic infections were caused by strains belonging to only 7 bacterial classes and one yeast genus. With less than 1% incidence in the majority of the countries, the taxonomic diversity of rare pathogens markedly increased adding another 13 bacterial classes to the list of aetiologies of sepsis. Even more, a survey of 54 case studies uncovered further unusual species that belonged to 27 classes. In order to break through the obstacles of poor cultural diagnosis, the use of culture-independent, pan-bacterial 16S rDNA PCR and sequencing analysis has entered routine as a helpful method of precise strain identification within hours instead of days.

Introduction

Latest 2016 recommendations of the Surviving Sepsis Campaign guide to the immediate (within one hour) administration of empiric broad-spectrum antimicrobials after diagnosis of sepsis or septic shock. The problem with gold standard blood culture is that cultures are often negative with blood draws after initiation of antibiotic treatment. Also, as demonstrated in a recent study (1), cultures grow positive with bacteria and yeasts only after 12 to 61 hours which delays targeted antibiotic treatment and increases the risk of poor prognosis of the patients. According to the recommendations, empiric antimicrobials should include combinations of two antimicrobials of different antibiotic

classes targeting the most probable, i.e., common pathogens.

Common pathogens

In clinical practice in the Western and Eastern hemisphere, *Escherichia coli*, *Staphylococcus aureus* and coagulase-negative staphylococci have been shown in clinical studies to make up 43.1 to 65% of cases of bacteraemia in sepsis (Table 1). In large studies from Germany (4380 patients), France (842), USA (1585), Saudi Arabia (1626) and Pakistan (521), further common pathogens or representatives of the genus were identified with incidences of sepsis cases between 1 and 10%, including streptococci, *Salmonella typhi/paratyphi*, *Pseudomonas* spp., *Acinetobacter* spp., *Klebsiella* spp., enterococci, *Candida* spp., *Enterobacter* spp., *Proteus* spp., *Serratia* spp., *Bacteroides* spp., *Clostridium* spp., and, in USA, *Mycobacterium* spp. (Table 1). Depending on the country, these 16 pathogens/groups cause sepsis in 32.9 to 47.6% of cases. Together, common pathogens account for 90.7 to 98.2% of systemic infections.

Rare pathogens

Other pathogenic bacteria were associated with infections at lower than 1% incidence. This arbitrarily defined group of rare pathogens accounted for 1.4 to 13.2% of bacteraemias in the studies of the 5 countries (Table 1). The pathogens represent 13 taxonomic classes, including Enterobacterales (*Citrobacter* spp., *Yersinia enterocolitica*), Pasteurellales (*Haemophilus influenzae*), Vibrionales (*Vibrio* spp.), Neisseriales (*Neisseria meningitidis*), Clostridiales (*Peptostreptococcus* spp., *Peptococcus* spp.), Rhizobiales (*Brucella* spp.), Burkholderiales (*Burkholderia* spp.), Propionibacterales (*Propionibacterium* spp.), Bacillales (*Listeria monocytogenes*, *Bacillus* spp.), Fusobacterales (*Fusobacterium* spp.), Campylobacterales (*Campylobacter fetus*), Corynebacterales (*Corynebacterium* spp.) and Veillonellales (*Veillonella* spp.). For comparison, the 16 common patho-

Table 1: Pathogens identified in blood cultures of septic patients.

| Strain | Culture incidence (%) | | | | | | Range | | Incidence in majority of countries |
|---|-----------------------|--------------|-------------|-------------|--------------|-------------|-------------|-------------|---------------------------------------|
| | Germany | France | USA | Pakistan | Saudi Arabia | lower | | | |
| | (2) | (3) | (4) | (5) | (6) | | | | |
| Coagulase-negative staphylococci | 9.3 | 10.7 | 9.2 | 36.4 | 24.7 | 9.2 | 36.4 | ≥10% | |
| <i>E. coli</i> | 23.9 | 29.2 | 15.0 | 18.9 | 9.1 | 9.1 | 29.2 | | |
| <i>Staphylococcus aureus</i> | 20.4 | 21.4 | 18.9 | 9.7 | 18.6 | 9.7 | 21.4 | | |
| Sum | 53.6 | 61.3 | 43.1 | 65.0 | 52.4 | 43.1 | 65.0 | | |
| <i>Streptococcus pneumoniae</i> and other <i>Streptococcus</i> spp. | 12.4 | 16.7 | 8.7 | 1.5 | 7.7 | 1.5 | 16.7 | ≥1% to <10% | |
| <i>Salmonella typhi/paratyphi</i> | 2.1 | 1.1 | | 8.9 | 0.7 | 0.7 | 8.9 | | |
| <i>Pseudomonas aeruginosa</i> and other <i>Pseudomonas</i> spp. | 4.9 | 3.8 | 6.3 | 4.4 | 8.7 | 4.4 | 8.7 | | |
| <i>Acinetobacter baumannii</i> and other <i>Acinetobacter</i> spp. | 2.2 | 1.1 | 1.4 | 3.4 | 7.8 | 1.1 | 7.8 | | |
| <i>Klebsiella pneumoniae</i> and other <i>Klebsiella</i> spp. | 5.4 | 4.3 | 6.9 | 4.0 | 6.7 | 4.0 | 6.9 | | |
| <i>Enterococcus faecalis</i> and other <i>Enterococcus</i> spp. | 5.5 | 5.5 | 6.7 | 5.7 | 4.8 | 4.8 | 6.7 | | |
| <i>Candida albicans</i> and other <i>Candida</i> spp. | | | 5.9 | 5.0 | 5.9 | 5.0 | 5.9 | | |
| <i>Enterobacter cloacae</i> and other <i>Enterobacter</i> spp. | 3.6 | | 2.6 | | 1.2 | 1.2 | 3.6 | | |
| <i>Proteus mirabilis</i> and other <i>Proteus</i> spp. | 2.3 | 2.7 | 1.7 | | 0.7 | 0.7 | 2.7 | | |
| <i>Serratia marcescens</i> and other <i>Serratia</i> spp. | 0.8 | | 2.3 | | 1.4 | 0.8 | 2.3 | | |
| <i>Bacteroides</i> spp. | 1.2 | | 1.7 | | 0.2 | 0.2 | 1.7 | | |
| <i>Clostridium perfringens</i> and other <i>Clostridium</i> spp. | 0.5 | | 1.6 | | | 0.5 | 1.6 | | |
| <i>Mycobacterium</i> spp. | | | 1.8 | | | | 1.8 | | |
| Sum | 40.9 | 35.2 | 47.6 | 32.9 | 45.8 | 32.2 | 47.6 | | |
| <i>Cryptococcus neoformans</i> | | | 0.8 | | | | 0.8 | <1% | |
| <i>Citrobacter</i> spp. | 0.7 | | | 0.6 | 0.3 | 0.3 | 0.7 | | |
| <i>Lactobacillus</i> spp. | 0.1 | | 0.6 | | | 0.1 | 0.6 | | |
| <i>Haemophilus influenzae</i> | 0.5 | | 0.3 | | 0.4 | 0.3 | 0.5 | | |
| <i>Stenotrophomonas maltophilia</i> | | | 0.5 | 0.2 | | 0.2 | 0.5 | | |
| <i>Vibrio</i> spp. | | | | 0.4 | | | 0.4 | | |
| <i>Neisseria meningitidis</i> | 0.3 | 0.2 | | | | 0.2 | 0.3 | | |
| <i>Peptostreptococcus</i> spp. | 0.3 | | | | | | 0.3 | | |
| <i>Yersinia enterocolitica</i> | 0.1 | | | 0.2 | | 0.1 | 0.2 | | |
| <i>Brucella</i> spp. | | | | | 0.2 | | 0.2 | | |
| <i>Burkholderia</i> spp. | | | | | 0.2 | | 0.2 | | |
| <i>Propionibacterium</i> spp. | 0.2 | | | | | | 0.2 | | |
| <i>Listeria monocytogenes</i> | 0.2 | | 0.1 | | | 0.1 | 0.2 | | |
| <i>Bacillus</i> spp. | 0.1 | | 0.1 | | 0.2 | 0.1 | 0.2 | | |
| <i>Fusobacterium</i> spp. | 0.1 | | | | | | 0.1 | | |
| <i>Campylobacter fetus</i> | 0.1 | | | | | | 0.1 | | |
| <i>Corynebacterium</i> spp. | 0.3 | | 0.1 | | | 0.1 | 0.3 | | |
| <i>Peptococcus</i> spp. | 0.1 | | | | | | 0.1 | | |
| <i>Veillonella</i> spp. | 0.02 | | | | | | 0.02 | | |
| <i>Aspergillus fumigatus</i> | 0.02 | | | | | | 0.02 | | |
| Other Gram-negatives (cumulated, not further specified) | 1.5 | 6.2 | 4.3 | | | 1.5 | 6.2 | | |
| Other Gram-positives (cumulated, not further specified) | | 1.9 | 0.6 | | 0.5 | 0.5 | 1.9 | | |
| Unspecified bacteria (cumulated) | | 3.2 | | | | | 3.2 | | |
| Unspecified fungi (cumulated) | 1.3 | 1.7 | | | | 1.3 | 1.7 | | |
| Sum | 5.7 | 13.2 | 7.6 | 1.4 | 1.8 | 1.4 | 13.2 | | |
| Total | 100.2 | 109.7 | 98.3 | 99.2 | 100.0 | | | | |

gens (Table 1) fall into only 7 classes (Bacillales, Enterobacteriales, Lactobacillales, Pseudomonadales, Bacteroidales, Clostridiales and Corynebacteriales).

There is a long record of further bacteria involved in rare cases of sepsis. A Google Scholar search under the term "rare pathogens in sepsis" resulted in 54 case reports uncovering an enormous diversity of unusual aetiologies of bacteraemia. Most of the identified bacteria normally live in the environment, on plants and in animals, and can give rise to sepsis in patients with neutropenia, meningitis, endocarditis, diabetic foot ulcers or catheterised, in the course of transplantations and surgeries, in neonates, after cat and dog bites and in other incidences. Table 2 sorts the pathogens into 40 families and 27 classes to underline the enormous diversity of bacterial species involved in septicaemia. This list does not claim to be comprehensive. Nonetheless, these case

reports emphasise the notion that unusual bacteria stem from taxonomically highly diverse groups.

This diversity of pathogenic bacteria can pose problems to the cultivation and the identification of strains in clinical routine, because standard media not always meet their growth conditions and phenotypic tests may not represent their characters for identification purposes. In fact, the diagnostic procedure of unexpected aetiologies of sepsis can be laborious and time-consuming and standard phenotypic identification practices may fail or lead to misidentification. For instance, in their report Malkan et al. (15) described a case of an anaemic patient supported with iron by a venous infusion port catheter. The patient developed sepsis by a methicillin-sensitive *S. aureus* which was treated with daptomycin. At day 3, the catheter was removed. Cultures of the device were positive with Gram-positive cocci and Gram-negative

Table 2: Cultured rare strains causing sepsis and detectability by 16S rRNA gene PCR and sequencing.

| Class / family | Species | Diagnosis | Reference | Strains found by direct 16S rDNA on taxonomic level ^a | No. species found (no. genera) ^b | Primers binding ^c |
|---------------------------|---|--|-----------|--|---|------------------------------|
| Actinomycetales | | | | | | |
| Actinomycetaceae | <i>Trueperella bernardiae</i> | sepsis, diabetic foot ulcer | (7) | Actinomycetaceae | 16 (3) | yes |
| Aeromonadales | | | | | | |
| Aeromonadaceae | <i>Aeromonas hydrophila</i> | paediatric sepsis | (8) | <i>A. hydrophila</i> | 1 | yes |
| Succinivibrionaceae | <i>Anaerobiospirillum succiniciproducens</i> | sepsis after cat bite | (9) | Aeromonadales | | yes |
| Alteromonadales | | | | | | |
| Shewanellaceae | <i>Shewanella algae</i> | neonatal sepsis | (10) | <i>Shewanella</i> | 2 | yes |
| Bacillales | | | | | | |
| Paenibacillaceae | <i>Paenibacillus amylolyticus</i> | bacteraemia | (11) | <i>Paenibacillus</i> | 2 | yes |
| Bacillaceae | <i>Lysinibacillus fusiformis/sphaericus</i> | bacteraemia | (11) | Bacillaceae | 79 (10) | yes |
| Bacteroidales | | | | | | |
| Prevotellaceae | <i>Prevotella buccae</i> | sepsis | (12) | <i>Prevotella</i> | 7 | yes |
| Burkholderiales | | | | | | |
| Burkholderiaceae | <i>Burkholderia cepacia</i> | neonatal sepsis | (13) | <i>B. cepacia</i> | 10 | yes |
| | <i>Ralstonia picketti</i> | neonatal sepsis | (14) | <i>R. picketti</i> | 5 | yes |
| Comamonadaceae | <i>Acidovorax avenae</i> | catheter-related | (15) | <i>A. avenae</i> | 6 | yes |
| | <i>Comamonas terrigena</i> | sepsis, endocarditis | (16) | <i>C. terrigena</i> | 5 | yes |
| | <i>Delftia acidovorans</i> | sepsis | (17) | <i>D. acidovorans</i> | 3 | yes |
| Bifidobacteriales | | | | | | |
| Bifidobacteriaceae | <i>Gardnerella vaginalis</i> | sepsis, post partum | (18) | <i>G. vaginalis</i> | 1 | yes |
| Campylobacteriales | | | | | | |
| Campylobacteraceae | <i>Campylobacter jejuni</i> | maternal sepsis | (19) | <i>C. jejuni</i> | 4 | yes |
| Caulobacteriales | | | | | | |
| Caulobacteraceae | <i>Brevundimonas vesicularis</i> | neonatal sepsis | (20) | <i>B. vesicularis</i> | 6 | yes |
| Chlamydiales | | | | | | |
| Chlamydiaceae | <i>Chlamydomphila abortus</i> | septicaemia | (21) | none | 0 | no result |
| Coriobacteriales | | | | | | |
| Atopobiaceae | <i>Atopobium</i> spp. | sepsis | (22) | <i>Atopobium</i> | 2 | yes |
| Corynebacteriales | | | | | | |
| Nocardiaceae | <i>Nocardia</i> spp. | sepsis | (23) | <i>Nocardia</i> | 11 | yes |
| | <i>Rhodococcus equi</i> | neonatal sepsis | (24) | <i>R. equi</i> | 4 | yes |
| Tsukamurellaceae | <i>Tsukamurella paurometabolum</i> | catheter-related | (25) | <i>Tsukamurella</i> | 6 | yes |
| Enterobacteriales | | | | | | |
| Enterobacteriaceae | <i>Cronobacter</i> spp. | neonatal sepsis | (26) | <i>Cronobacter</i> | 2 | yes |
| | <i>Leclercia adecarboxylata</i> | sepsis | (27) | Enterobacteriaceae | 71 (19) | yes |
| | <i>Shigella sonnei</i> | sepsis, neutropaenia | (28) | <i>S. sonnei</i> | 4 | yes |
| Hafniaceae | <i>Edwardsiella tarda</i> | neonatal sepsis | (29) | Hafniaceae | 1 | yes |
| | <i>Hafnia alvei</i> | sepsis | (30) | <i>H. alvei</i> | 1 | yes |
| Morganellaceae | <i>Morganella morganii</i> | neonatal sepsis | (31) | <i>M. morganii</i> | 1 | yes |
| | <i>Providencia stuartii</i> | sepsis | (32) | <i>P. stuartii</i> | 3 | yes |
| Erwiniaceae | <i>Pantoea dispersae</i> | neonatal sepsis | (33) | <i>P. dispersae</i> | 6 | yes |
| Yersiniaceae | <i>Rahnella aquatilis</i> | sepsis, contaminated intravenous fluid | (34) | <i>R. aquatilis</i> | 1 | yes |
| Unclassified | <i>Plesiomonas shigelloides</i> | sepsis, meningitis | (35) | <i>P. shigelloides</i> | 1 | yes |
| Erysipelotrichales | | | | | | |
| Erysipelotrichaceae | <i>Erysipelothrix rhusiopathiae</i> | sepsis, endocarditis | (36) | Erysipelotrichaceae | 1 | no result |
| Flavobacteriales | | | | | | |
| Flavobacteriaceae | <i>Capnocytophaga canimorsus</i> | sepsis, dog bite | (37) | <i>Capnocytophaga</i> | 2 | yes |
| | <i>Chryseobacterium indologenes</i> | neonatal sepsis, meningitis | (38) | <i>C. indologenes</i> | 5 | yes |
| | <i>Elizabethkingia (Chryseobacterium) meningosepticum</i> | sepsis, diabetic patient | (39) | <i>E. meningosepticum</i> | 1 | yes |
| | <i>Empedobacter brevis</i> | neonatal sepsis | (40) | Flavobacteriaceae | 12 (6) | yes |
| Fusobacteriales | | | | | | |
| Fusobacteriaceae | <i>Fusobacterium necrophorum</i> | sepsis | (41) | <i>F. necrophorum</i> | 7 | yes |
| Leptotrichiaceae | <i>Leptotrichia trevisanii</i> | sepsis, | (42) | <i>Leptotrichia</i> | 1 | yes |
| Lactobacillales | | | | | | |
| Leuconostocaceae | <i>Leuconostoc</i> spp. | neonatal sepsis | (43) | <i>Leuconostoc</i> | 2 | no result |
| | <i>Weissella confusa</i> | bacteraemia | (44) | <i>W. confusa</i> | 4 | yes |
| Micrococcales | | | | | | |
| Micrococaceae | <i>Kocuria kristinae</i> | paediatric sepsis | (45) | <i>K. kristinae</i> | 7 | yes |
| | <i>Rothia dentocariosa</i> | paediatric sepsis | (46) | <i>R. dentocariosa</i> | 5 | yes |
| Promicromonosporaceae | <i>Cellulosimicrobium cellulans</i> | neonatal sepsis | (47) | <i>Cellulosimicrobium</i> | 1 | yes |
| Neisseriales | | | | | | |
| Chromobacteriaceae | <i>Chromobacterium violaceum</i> | paediatric sepsis | (48) | none | 0 | yes |
| Neisseriaceae | <i>Eikenella corrodens</i> | postanginal sepsis | (49) | Neisseriaceae | 19 (4) | yes |
| Pasteurellales | | | | | | |
| Pasteurellaceae | <i>Pasteurella multocida</i> | neonatal sepsis, meningitis | (50) | <i>P. multocida</i> | 5 | yes |

Table 2: Continued.

| Class / family | Species | Diagnosis | Reference | Direct 16S rDNA result on taxonomic level ^a | No. species found (no. genera) ^b | Primers binding ^c |
|----------------------------|--|-----------------------------|-----------|--|---|------------------------------|
| Pseudomonadales | | | | | | |
| Moraxellaceae | <i>Moraxella catarrhalis</i> | sepsis, cellulitis | (51) | <i>M. catarrhalis</i> | 5 | yes |
| Rhizobiales | | | | | | |
| Rhizobiaceae | <i>Agrobacterium tumefaciens (radiobacter)</i> | catheter-related sepsis | (52) | <i>A. tumefaciens</i> | 2 | yes |
| Brucellaceae | <i>Ochrobactrum anthropi</i> | biliary sepsis | (53) | <i>O. anthropi</i> | 7 | yes |
| Sphingomonadales | | | | | | |
| Sphingomonadaceae | <i>Sphingomonas paucimobilis</i> | paediatric sepsis | (54) | <i>S. paucimobilis</i> | 22 | yes |
| Tissierellales | | | | | | |
| Peptoniphilaceae | <i>Peptoniphilus lacrimalis</i> | sepsis, diabetic foot ulcer | (7) | <i>P. lacrimalis</i> | 10 | yes |
| Veillonellales | | | | | | |
| Veillonellaceae | <i>Veillonella parvula</i> | sepsis | (55) | <i>V. parvula</i> | 5 | no result |
| Vibrionales | | | | | | |
| Vibrionaceae | <i>Vibrio vulnificus</i> | septicaemia | (56) | <i>V. vulnificus</i> | 2 | yes |
| Unclassified | | | | | | |
| Gammaproteobacteria | <i>Ignatzschineria indica</i> | sepsis | (57) | none | 0 | yes |
| | <i>Wohlfahrtiimonas chitiniclastica</i> | sepsis | (57) | none | 0 | yes |

^a Results from evaluations of Molzym's products involving other diseases and clinical fluid and tissue materials; products tested: SepsitTest™-UMD CE, UMD-SelectNA™ CE or Micro-Dx™ CE (all include DNA extraction and PCR assay); amplicons from positive tests were Sanger-sequenced and Blast-analysed (sepsitest-blast.net and https://blast.ncbi.nlm.nih.gov/Blast.cgi?PAGE_TYPE=BlastSearch)

^b See also the complete list of organisms detected directly in clinical and other materials (https://www.molzym.com/images/products/Flyer_App_Notes/molzym_list-of-microbes_03_2018.pdf)

^c Binding to 16S rRNA gene of cultured strain; Primer designing tool, NCBI (<https://www.ncbi.nlm.nih.gov/tools/primer-blast/>); specificity stringency: target sequences with 6 or more mismatches to forward or reverse primers or with 2 or more mismatches at 3' end were excluded (no result); 16S rRNA gene sequences were retrieved from NCBI (<https://www.ncbi.nlm.nih.gov/gene/>)

rods. Same day blood cultures grew only Gram-negative rods and treatment was changed to gentamicin. Subculturing was performed on various media resulting in very small colonies, and Vitek 2 (bioMérieux) compact Gram-negative identification resulted in *Cupriavidus pauculus*. The bioMérieux API 20E system could not assign the organism to a certain species. At the end, definitive identification of the strain as *Acidovorax avenae* was reached by sequencing analysis of a 485 bp PCR amplicon of the 16S rRNA gene.

Ribosomal RNA gene sequencing in clinical practice

Ribosomal RNA gene sequencing is used in routine laboratories as a reliable method of strain identification. Mahlen and Clarridge (58), for instance, elaborated a strategy of broad-range 16S rRNA gene sequencing analysis of Gram-negative rods and coccobacilli as aid for routine identification methods. For strain identification, the rRNA genes or parts thereof are amplified in a PCR reaction. The reaction includes PCR primers that bind to sequences highly conserved among the kingdom of bacteria, while the sequences in-between vary among species. Identification of an unknown strain is achieved by comparison of its determined sequence to a library of sequences of known strains. Online tools are available that perform the assignment by an elaborated algorithm. Examples are the NCBI Basic Local Alignment Search Tool

(BLAST®) leBIBI and Molzym's SepsitTest™BLAST.

Culture-independent diagnosis by rRNA gene sequencing

rRNA gene sequencing approaches, including Sanger and NGS sequencing of amplified and cloned 16S rRNA gene sequences, respectively, have been used for the direct, i.e., culture-independent analysis of bacterial strains and multiple infections for quite a while. Many of the reported methods are used for research purposes to answer specific scientific questions rather than for use of surveillance and in daily routine diagnostic practice.

rRNA gene sequencing promises to add benefit to culture diagnosis in two clinical situations: diagnosis of i) fastidious and anaerobic bacteria with unusual growth requirements and ii) growth-inhibited, but not cleared pathogens as a result of antimicrobial treatment of patients. Reported applications of the use of rRNA gene sequencing comprise of diseases like pneumonia, meningitis, thrombosis, pleuritis, septic arthritis, septic coxitis, intra-abdominal infection, intra-amniotic infection, brain abscess, recurrent cholecystitis, spondylitis, osteomyelitis, septicaemia, sepsis, subcutaneous abscess, spinal abscess, empyema (69), endophthalmitis (70), cystic fibrosis (71), infective endocarditis (72), urinary tract infections (73), airway inflammation (74), chronic venous leg ulcer (75), prosthetic joint infections (76) and others. Unlike other as-

says for direct analysis of a limited number of selected common pathogens, the unbiased nature of broad-range 16S rRNA gene PCR and sequencing follows a universal approach of the identification of common, rare and unusual pathogens.

Standardisation

The majority of methods for direct rRNA gene sequencing analysis of pathogens in clinical specimens is in-house developed tests and thus are not comparable among laboratories. Such tests differ in their pre-analytical and analytical procedures. A variety of commercial kits used in the laboratories are usually for DNA extraction from clinical samples and amplification. The problem is that chemicals, reagents and consumables are generally not produced under control of the absence of microbial DNA (77), because they are designed for applications where DNA contamination is irrelevant. Under these conditions the avoidance of false positive results in the diagnosis of low bacterial load materials is difficult to reach if not impossible. Moreover, the sensitivity and specificity of tests can be lowered by the vast load of human DNA in many clinical specimens (78). Another consequence is that Sanger and NGS sequencing analysis can be disturbed by human sequences from unspecific amplification (78).

The only unbiased 16S rRNA gene PCR and sequencing tests, which are CE-approved for the in-vitro diagnostics of hundreds of pathogens are provided by Molzym. The products are complete solutions for the whole process, including enrichment (depletion of human DNA) and extraction of DNA from bacteria present in samples at low loads (<100 cfu/ml), detection by Real-Time PCR amplification and sequencing analysis of the amplified hypervariable V3/V4 region of the 16S rRNA gene. By employing highly active and contamination-free PCR reagents, amplification runs of 40 cycles are performed without signals in negative PCR controls. Solutions are available for manual extraction, SepsiTest™-UMD, semi-automated extraction, UMD-SelectNA™, and fully automated extraction, Micro-Dx™ in the SelectNA™*plus* robot. Other, contamination-free research-use-only products provide options for the extraction of bacterial DNA from fluid (MolYsis™ kits) and tissue samples (Ultra-Deep Microbiome kit) and amplification of custom targets as well as defined parts of the 16S rRNA gene for qualitative, quantitative Real-Time PCR (Mastermix 16S/18S kits) and library preparation purposes (NGSeq 16S V3/V4). Validations in clinical studies and routine use have proven that unbiased 16S rRNA gene analysis can reliably detect and identify a great variety of microorganisms, in-

cluding rare pathogens directly in fluid and tissue biopsies (59-68).

All of the common and rare species found in the studies discussed above (Table 1) have been shown to be detectable in other evaluations. Among the 54 unusual species listed in Table 2, 31 have been identified by Molzym's direct 16S rRNA gene sequencing approaches in studies employing other specimens. In 19 of the cases, 16S rRNA gene sequencing evaluations have not yet identified the same strains, but other species of the same genus, family or class of the cultured strains. No relatives of *Chlamydophila abortus*, *Chromobacterium violaceum*, *Ignatzschineria indica* and *Wohlfahrtimonas chitinoclastica* have been found at the family or class level so far (Table 2). However, except for *C. abortus*, *in silico* testing of primer binding to the 16S rRNA gene targets was positive, indicating that the strains are detectable by the assay. Thus, 16S rRNA gene sequencing as represented by Molzym's molecular diagnostic kits have the potential to identify essentially all pathogens directly found in samples. The kits are CE-marked for the *in-vitro* diagnosis of pathogens.

Conclusions

Ribosomal RNA gene amplification and sequencing analysis is a frequently used method for the unbiased identification of rare and unusual pathogenic bacteria in clinical routine. In the past years, the method has been adapted to the culture-independent diagnosis of infected clinical materials at low loads of pathogens. With this approach, however, a variety of problems have been faced that disturb the precise diagnosis. Factors negatively influencing sensitivity, specificity and sequencing identification include, among others contamination of extraction and amplification reagents and consumables with microbial DNA and a great excess of human DNA. Molzym addressed these problems by the supply of contamination-free kits and reagents that enable the processing of human DNA-depleted preparations from a variety of clinical samples by just one protocol and the amplification over 40 cycles without background in negative PCR controls. By this solution hundreds of rare and unusual pathogens can be identified within hours without the need of cultivation.

References

1. **Ning Y, Hu R, Yao G, Bo S (2016)** Time to positivity of blood culture and its prognostic value in bloodstream infection. *Eur J Clin Microbiol Inf Dis* 35, 619–624.

2. **Rosenthal EJK (1993)** Epidemiologie von Septikämie-Erregern. Dtsch Med Wschr 118, 1269-1275 (in German).
3. **Brun-Buisson C, Doyon F, Carlet J et al. (1996)** Bacteremia and severe sepsis in adults: a multicentre prospective survey in ICUs and wards of 24 hospitals. Am J Respir Crit Care Med 154, 617-624.
4. **Weinstein MP, Towns ML, Quartey SM, Mirrett S, Reimer LG, Parmigiani G, Relier LB (1997)** The clinical significance of positive blood cultures in the 1990s: a prospective comprehensive evaluation of the microbiology, epidemiology, and outcome of bacteremia and fungemia in adults. Clin Inf Dis 24, 584-602.
5. **Sohail M, Sultana Q, Rasool K, Sarwar S, Basit A, Khalid M (2015)** Bacteremia prediction by inflammatory factors and recent trend in drug resistance of bacteria isolated from blood stream infection. J Inf Mol Biol 3, 75-80.
6. **Asghar AH (2006)** Frequency and antimicrobial susceptibility patterns of bacterial pathogens isolated from septicemic patients in Makkah hospitals. Saudi Med J 27, 361-367.
7. **Schneider U, Ekenberg C, Sode N, Knudsen J (2015)** A case of diabetic foot ulcers complicated by severe infection and sepsis with *Trueperella bernardiae*. JMM Case Rep 2, doi:10.1099/jmmcr.0.000006.
8. **Rodríguez CN, Campos R, Pastran B, Jimenez I, Garcia A, Meijomil P, Rodríguez-Morales AJ (2005)** Sepsis due to extended-spectrum β -lactamase-producing *Aeromonas hydrophila* in a pediatric patient with diarrhea and pneumonia. Clin Inf Dis 41, 421-422.
9. **I.M. Markusse IM, van Mens SP, de Vries PJ (2018)** An unusual causative pathogen of sepsis after a cat bite: *Anaerobiospirillum succiniciproducens*. Netherlands J Med 76, 87-89.
10. **Charles MV, Srirangaraj S, Kali A (2015)** Neonatal sepsis caused by *Shewanella algae*: a case report. Australas Med J 8, 64-6.
11. **Wenzler E, Kamboj K, Balada-Llasat JM (2015)** Severe sepsis secondary to persistent *Lysinibacillus sphaericus*, *Lysinibacillus fusiformis* and *Paenibacillus amylolyticus* bacteremia. Intern J Inf Dis 35, 93-95.
12. **Jensen KS, Biggs KA, Cardwell MS (2015)** Retropharyngeal abscess complicated by *Prevotella buccae* sepsis during pregnancy. A case report. J Reprod Med 60, 87-89.
13. **Patra S, Bhat YR, Lewis LE, Purakayastha J, Sivaramaraju VV, Kalwaje E V, Mishra S (2014)** *Burkholderia cepacia* sepsis among neonates. Indian J Pediatr 81, 1233, <https://doi.org/10.1007/s12098-014-1473-9>.
14. **Sharma D, Sharma P, Soni P, Gupta B (2017)** *Ralstonia picketti* neonatal sepsis: a case report. BMC Res Notes 10, 28, <https://doi.org/10.1186/s13104-016-2347-1>.
15. **Malkan AD, Strollo W, Scholand SJ, Dudrick SJ (2009)** Implanted-port-catheter-related sepsis caused by *Acidovorax avenae* and methicillin-sensitive *Staphylococcus aureus*. J Clin Microbiol 47, 358-3361.
16. **Sonnenwirth AC (1970)** Bacteremia with and without meningitis due to *Yersinia enterocolitica*, *Edwardsiella tarda*, *Comamonas terrigena*, and *Pseudomonas maltophilia*. Ann New York Acad Sci 174, 488-502.
17. **Bilgin H, Sarmis A, Tigen E, Soyletir G, Mulazimoglu L (2015)** *Delftia acidovorans*: a rare pathogen in immunocompetent and immunocompromised patients. Can J Inf Dis Med Microbiol 26, 277-279.
18. **Murray L, Halpin J, Casserly B, O'Connell NH, Scanlon T (2019)** A pyohydropneumothorax with sepsis, secondary to *Gardnerella vaginalis* infection in a postpartum female. Resp Med Case Rep 26, 189-192.
19. **Meyer A, Stallmach T, Goldenberger D, Altwegg M (1997)** Lethal maternal sepsis caused by *Campylobacter jejuni*: pathogen preserved in placenta and identified by molecular methods. Mod Pathol 10, 1253-1256.
20. **Karadag N, Karagol BS, Dursun A, Okumus N, Tanir G, Zenciroglu A (2012)** A premature neonate with early-onset neonatal sepsis owing to *Brevundimonas vesicularis* complicated by persistent meningitis and lymphadenopathy. Paed Inter Child Health 32, 239-241.
21. **Walder G, Hotzel H, Brezinka C, Gritsch W, Tauber R, Würzner R, Ploner F (2005)** An unusual cause of sepsis during pregnancy: recognizing infection with *Chlamydophila abortus*. Obst Gynecol 106, 1215-1217.
22. **Salimnia H, Noronha A, Sobel JD, Fairfax MR (2008)** Sepsis associated with a new

- Atopobium* species, provisionally named *Atopobium detroitii*: case report and review of the current status of the species *Atopobium*. Scand J Inf Dis 40, 679-681.
23. **Cassar C L (2007)** *Nocardia* sepsis in a multigravida with systemic lupus erythematosus and autoimmune hepatitis. Anaesth Int Care 35, 601-614.
 24. **Bildik HN, Takcı S, Yurdakök MKara A (2013)** Neonatal sepsis due to *Rhodococcus equi* in two preterm infants. Turk J Ped 55, 229-231.
 25. **Lai KK (1993)** A cancer patient with central venous catheter-related sepsis caused by *Tsukamurella paurometabolum* (*Gordona aurantiaca*). Clin Inf Dis 17, 285-287.
 26. **Hunter CJ, Bean JF (2013)** *Cronobacter*: an emerging opportunistic pathogen associated with neonatal meningitis, sepsis and necrotizing enterocolitis. J Perinatol 33, 581-585.
 27. **Matsuura H, Sugiyama S (2018)** Sepsis and *Leclercia adecarboxylata*. QJM: Intern J Med 111,733-734.
 28. **Kenet G, Salomon F, Samra Z, Pinkhas J, Sidi Y, Arber N (1994)** Fatal *Shigella* sepsis in a neutropenic patient. Mt Sinai J Med 61, 367-368.
 29. **Hashavya S, Averbuch D, Berger I, Ofek-Shlomai N, Pitashny M, Hidalgo C, Ergaz Z (2011)** Neonatal sepsis following maternal amnionitis by *Edwardsiella tarda*: a case report and a review of the literature. Eur J Pediatr 170: 111-113.
 30. **Casanova-Román M, Sanchez-Porto A, Casanova-Bellido M (2004)** Late-onset neonatal sepsis due to *Hafnia alvei*. Scan J Inf Dis 36, 70-71.
 31. **Salen PN, Eppes S (1997)** *Morganella morganii*: a newly reported, rare cause of neonatal sepsis. Acad Emerg Med 4, 711-714.
 32. **Avolio M, Diamante P, Zamparo S, Modolo ML, Grosso S, Zigante P, Tosoni N, De Rosa R, Stano P, Camporese A (2010)** Molecular identification of bloodstream pathogens in patients presenting to the emergency department with suspected sepsis. Shock 34, 27-30.
 33. **Veerendra M, Dinesh Y, Jyoti S, Nidhi G, Kuldeep S (2013)** *Pantoea dispersa*: an unusual cause of neonatal sepsis. Braz J Infect Dis 17, 726-728.
 34. **Chang CL, Jeong J, Shin JH, Lee EY, Son HC (1999)** *Rahnella aquatilis* sepsis in an immunocompetent adult. J Clin Microbiol 37, 4161-4162.
 35. **Stock I (2004)** *Plesiomonas shigelloides*: an emerging pathogen with unusual properties. Rev Med Microbiol 15, 129-139.
 36. **Reboli AC, Farrar WE (1989)** *Erysipelothrix rhusiopathiae*: an occupational pathogen. Clin Microbiol Rev 2, 354-359.
 37. **Butler T (2015)** *Capnocytophaga canimorsus*: an emerging cause of sepsis, meningitis, and post-splenectomy infection after dog bites. Eur J Clin Microbiol Infect Dis 34, 1271-1280.
 38. **Eshwara VK, Sasi A, Munim F, Purkayastha J, Lewis LE, Mukhopadhyay C (2014)** Neonatal meningitis and sepsis by *Chryseobacterium indologenes*: a rare and resistant bacterium. Indian J Pediatr 81, 611-613.
 39. **Lee SW, Tsai CA, Lee BJ (2008)** *Chryseobacterium meningosepticum* sepsis complicated with retroperitoneal hematoma and pleural effusion in a diabetic patient. J Chin Med Assoc 71, 473-476.
 40. **Basani L, Aepala R (2018)** *Empedobacter brevis* causing early onset sepsis and pneumonia in a neonate: case report and review of literature. Int J Contemp Ped 5, 654-656.
 41. **Tawa A, Larmet R, Malledant Y, Seguin P (2016)** Severe sepsis associated with Lemierre's syndrome: a rare but life-threatening disease. Case Rep Crit Care, <https://doi.org/10.1155/2016/1264283>.
 42. **Schrimsher JM, McGuirk JP, Hinthorn DR (2013)** *Leptotrichia trevisanii* sepsis after bone marrow transplantation. Emerg Infect Dis 19, 1690-1691.
 43. **Yossuck P, Miller-Canfield P, Moffett K, Graeber J (2009)** *Leuconostoc* spp sepsis in an extremely low birth weight infant: a case report and review of the literature. West Virginia Med J 105, <https://link.galegroup.com/apps/doc/A217847901/AONE?u=googlescholar&sid=AONE&xid=f89def74>.
 44. **Kumar A, Augustine D, Sudhindran S, Kurian AM, Dinesh KR, Karim S, Philip R (2015)** *Weissella confusa*: a rare cause of vancomycin-resistant Gram-positive bacteraemia. J Med Microbiol 60, 1539-1541.
 45. **Chen HM, Chi H, Chiu NC, Huang FY (2015)** *Kocuria kristinae*: a true pathogen in pediatric patients. J Microbiol Immunol Inf 48, 80-84.

46. Wiesmayr S, Stelzmueller I, Berger N, Jungraithmayr TC, Fille M, Eller M, Zimmerhackl LB, Margreiter R, Bonatti H (2006) *Rothia dentocariosa* sepsis in a pediatric renal transplant recipient having post-transplant lymphoproliferative disorders. *Ped Transplant* 10, 377-379.
47. Casanova-Román M, Sanchez-Porto A, Gomar JL, Casanova-Bellido M (2010) Early-onset neonatal sepsis due to *Cellulosimicrobium cellulans*. *Infection* 38, 321-323.
48. Vijayan AP, Anand MR, Remesh P (2009) *Chromobacterium violaceum* sepsis in an infant. *Indian Ped* 46, 721-722.
49. Celikel TH, Muthuswamy PP (1984) Septic pulmonary emboli secondary to internal jugular vein phlebitis (postanginal sepsis) caused by *Eikenella corrodens*. *Am Rev Res Dis* 130, <https://www.atsjournals.org/doi/abs/10.1164/arrd.1984.130.3.510>.
50. Spadafora R, Pomerio G, Delogu A, Gozzoli L, Gancia P (2011) A rare case of neonatal sepsis/meningitis caused by *Pasteurella multocida* complicated with status epilepticus and focal cerebritis. *Ped Med Chir* 33, 199-202.
51. Westendorp IC, Tiemessen MA, de Jong M, Soomers A, Wakelkamp IM, Boersma WG (2005) *Moraxella catarrhalis* sepsis in a patient with juvenile spinal muscle atrophy. *Neth J Med* 63,227-229.
52. Jeong SJ, Choi HK, Lee HW, Shin SH, Kim CK, Lee HS, Chin BS, Choi SH, Han SH, Kim MS, Kim CO, Choi JY, Song YG, Kim JM (2007) *Infect Chemother* 39, 270-273.
53. Wi YM, Peck KR (2019) Biliary sepsis caused by *Ochrobactrum anthropi*. *Jpn J Infect Dis* 63, 444-446.
54. Tat F, Devrim I, Gülfidan G, Gunay I, Dizdärer C (2012) A rare etiologic agent of sepsis in children: *Sphingomonas paucimobilis*. *Inf Dis Clin Practice* 20, 152–153.
55. Strach M, Siedlar M, Kowalczyk D, Zembala M, Grodzicki T (2006) Sepsis caused by *Veillonella parvula* infection in a 17-year-old patient with X-linked agammaglobulinemia (Bruton's disease). *J Clin Microbiol* 44, 2655-2656.
56. Shiang-Fu Huang SF, Chen CT (2004) *Vibrio vulnificus* — a rare but fulminant pathogen causing airway obstruction. *Am J Emerg Med* 23, 221–222.
57. Lysaght T, Wooster ME, Jenkins PC, Koniaris LG (2018) Myiasis-induced sepsis – a rare case report of *Wohlfahrtimonas chitinoclastica* and *Ignatzschineria indica* bacteremia in the continental United States. *Medicine* 97, p e13627, doi: 10.1097/MD.00000000000013627.
58. Mahlen SD, Clarridge JE (2011) Evaluation of a selection strategy before use of 16S rRNA gene sequencing for the identification of clinically significant Gram-negative rods and coccobacilli. *Am J Clin Pathol* 136, 381–388.
59. Stavnsbjerg C, Frimodt-Møller N, Moser C, Bjarnsholt T (2017) Comparison of two commercial broad range PCR and sequencing assays for identification of bacteria in culture-negative clinical samples. *BMC Inf Dis* 17, 233 DOI 10.1186/s12879-017-2333-9.
60. Peeters B, Herijgers P, Beuselinck K, Verhaegen J, Peetermans WE, Herregods MC, Desmet S, Lagrou K (2017) Added diagnostic value and impact on antimicrobial therapy of 16S rRNA PCR and amplicon sequencing on resected heart valves in infective endocarditis: a prospective cohort study. *Clin Microbiol Inf* 23, 888.e1–888.e5 (doi: <http://dx.doi.org/10.1016/j.cmi.2017.06.008>)
61. Nieman AE, Savelkoul PHM, Beishuizen A, Henrich B, Lamik B, MacKenzie CR, Kindgen-Milles D, Helmers A, Diaz C, Sakka SG Schade RP (2016) A prospective multicenter evaluation of direct molecular detection of blood stream infection from a clinical perspective. *BMC Inf Dis* 16, 314 DOI 10.1186/s12879-016-1646-4.
62. Trinick RE, Bunni L, Thorburn K, Hackett AP, Dalzell M, McNamara PS (2015) An observational study examining the relationship between respiratory symptoms, airway inflammation and bacteriology in children with severe neurodisability. *PLoS ONE* 10(4): e0124627. doi:10.1371/journal.pone.0124627.
63. Marsch G, Orszag P, Mashaqi B, Kuehn C, Haverich A (2015) Antibiotic therapy following polymerase chain reaction diagnosis of infective endocarditis: a single centre experience. *Interac Cardiovasc Thorac Surg* 20, 589-593.
64. Rogina P, Skvarc M, Stubljär D, Kofol R, Kaasch A (2014) Diagnostic utility of broad range bacterial 16S rRNA gene PCR with degradation of human and free bacterial

- DNA in bloodstream infection is more sensitive than an in-house developed PCR without degradation of human and free bacterial DNA. *Mediators of Inflammation*, doi.org/10.1155/2014/108592.
65. **Orszag P, Disqué C, Keim S, Lorenz MG, Wiesner O, Hadem J, Stiesch M, Haverich A, Kühn C (2014)** Monitoring of patients supported by extracorporeal membrane oxygenation for systemic infections by broad-range rRNA gene PCR amplification and sequence analysis. *J Clin Microbiol* 52, 307-311.
 66. **Haag H, Locher F, Nolte O (2013)** Molecular diagnosis of microbial aetiologies using SepsisTest™ in the daily routine of a diagnostic laboratory. *Diag Microbiol Inf Dis* 76, 413-418.
 67. **Grif K, Heller W., Proding M, Lechleitner K, Lass-Flörl C, and Orth D (2012)** Improvement of detection of bacterial pathogens in normally sterile body sites with a focus on orthopedic samples by use of a commercial 16S rRNA broad-range PCR and sequence analysis. *J Clin Microbiol* 50, 2250-2254.
 68. **Wellinghausen N, Kochem AJ, Disqué C, Mühl H, Gebert S, Winter J, Matten J, Sakka SG (2009)** Diagnosis of bacteremia in whole-blood samples by use of a commercial universal 16S rRNA gene-based PCR and sequence analysis. *J Clin Microbiol* 47, 2759-2765.
 69. **Rantakokko-Jalava K, Nikkari S, Jalava J, Eerola E, Skurnik M, Meurman O, Ruuskanen O, Alanen A, Kotilainen E, Toivanen P, Kotilainen P (2000)** Direct amplification of rRNA genes in diagnosis of bacterial infections. *J Clin Microbiol* 38, 32-39.
 70. **C.Michele Knox CM, Vicky Cevallos V, Todd P Margolis TP, Deborah Dean D (1999)** Identification of bacterial pathogens in patients with endophthalmitis by 16S ribosomal DNA typing. *Am J Ophthalmol* 128, 511-512.
 71. **Salipante SJ, Sengupta DJ, Rosenthal C, Costa G, Spangler J, et al. (2013)** Rapid 16S rRNA next-generation sequencing of polymicrobial clinical samples for diagnosis of complex bacterial infections. *PLOS ONE* 8(5): e65226. <https://doi.org/10.1371/journal.pone.0065226>.
 72. **Marsch G, Orszag P, Mashaqi B, Kuehn C, Haverich A (2015)** Antibiotic therapy following polymerase chain reaction diagnosis of infective endocarditis: a single centre experience. *Interac Cardiovasc Thorac Surg* 20, 589-593.
 73. **Sabat AJ, van Zanten E, Akkerboom V, Wisselink G, van Slochteren K, de Boer RF, Hendrix R, Friedrich AW, Rossen JWA, Kooistra-Smid AMD (2017)** Targeted next-generation sequencing of the 16S-23S rRNA region for culture-independent bacterial identification - increased discrimination of closely related species. *Scientific Reports* 7, article number: 3434.
 74. **Trinick RE, Bunni L, Thorburn K, Hackett AP, Dalzell M, McNamara PS (2015)** An observational study examining the relationship between respiratory symptoms, airway inflammation and bacteriology in children with severe neurodisability. *PLoS ONE* 10(4): e0124627. doi:10.1371/journal.pone.0124627.
 75. **Thomsen TR, Aasholm MS, Rudkjøbing VB, Saunders AM, Bjarnsholt T, Givskov M, Kirketerp-Møller K, Nielsen PH (2010)** The bacteriology of chronic venous leg ulcer examined by culture-independent molecular methods. *Wound Repair Regeneration* 18: 38-49.
 76. **Grif K, Heller W., Proding M, Lechleitner K, Lass-Flörl C, and Orth D (2012)** Improvement of detection of bacterial pathogens in normally sterile body sites with a focus on orthopedic samples by use of a commercial 16S rRNA broad-range PCR and sequence analysis. *J Clin Microbiol* 50, 2250-2254.
 77. **Lorenz MG (2019)** DNA-free reagents and materials for molecular pathogen analysis. *Res Mol Microbiol* 1/19: 1-9.
 78. **Handschr M, Karlic H, Hertl C, Pfeilstöcker M, Haslberger AG (2009)** Preanalytic removal of human DNA eliminates false signals in general 16S rDNA PCR monitoring of bacterial pathogens in blood. *Comp Immunol Microbiol Infect Dis* 32, 207-219.