



## Inactivation Type

(special for COVID test)



### Safe

Medical materials to prevent leakage and ensure biosafety



### Stable

BSA to form protective film of the virus to ensure its integrity



### Effective

Antibiotics are used to maintain the preservation of the specimen





### Convenient

Medical materials to prevent leakage and ensure biosafety


# Disposable Virus Sampling Kit

## Flocked swab advantages


- Non-toxic to cells and virus, RNase-free, enhance diagnostic sensitivity
- Complete and fast release of samples (up to 90%)




## Operation Procedures




1  
Use sampling swab to collect specimens




2  
Place the swab with collected specimen in the virus transport media.



3  
Break the swab at the breaking point, and discard the tail



4  
Tighten the lid of the media.



5  
Mark the information on the tube.

## Product Information

Product Code	Virus Transport Media Type	Tube Volume	Solution Volume	Swab Type
LVT01-1	Inactivation	5ml	2ml	Flocked oral swab
LVT01-2	Non-inactivation	5ml	2ml	Flocked oral swab
LVT02-1	Inactivation	10ml	3ml	Flocked oral swab
LVT02-2	Non-inactivation	10ml	3ml	Flocked oral swab
LVT03-1	Inactivation	10ml	3ml	Flocked nasal swab
LVT03-2	Non-inactivation	10ml	3ml	Flocked nasal swab
LVT04-1	Inactivation	10ml	3ml	Flocked oral & nasal swab
LVT04-2	Non-inactivation	10ml	3ml	Flocked oral & nasal swab



**ZHUHAI LITUO BIOTECHNOLOGY CO., LTD.**

Add : No.35, Yongan Three Road, Hongqi Town, Jinwan District, Zhuhai, Guangdong, China.  
Tel.: +86 756 8639521 (Sales Cenrter) FAX: +86 756 8639155 Email : market@lituo.com.cn



**HUNAN LITUO BIOTECHNOLOGY CO.,LTD.**

Add: No. 129, Weiluoke Road, Wangcheng Economic Development Zone, Changsha, Hunan, China.  
Tel.: +86 731 82788178 Fax: +86 731 82789603 Email: market@lituo.com.cn

