

MedSchenker® Vial Tray

Product Insert & How-to-Use Guide

INTENDED USE

MedSchenker® Vial Tray is used to hold and transport MedSchenker® Smart Transport Medium vials as well as other skirted vials for specimen transport and storage in quantities of up to 100 vials at a time.

FEATURES

- 219 x 182 x 20 mm rectangle shape with 15.6 x 20 mm wells
- Made of 100% polypropylene;
- Holds 100 vials
- Wells are positioned in a staggered formation to increase the ease in which vials are removed
- Trays have a low profile which helps save space
- Wells have an opening in the bottom to facilitate drainage
- Rectangle shaped tray has side cut-outs to make transport easy
- Trays can be stacked, when not filled
- Can be stored at temperatures between 2 and 25 degrees Celsius

ASSEMBLY

No assembly required

Manufacturer: Korea Standard
Jeollanamdo, South Korea

Exclusive Distributor:

Industrieweg 24C

4153 BW Beesd

The Netherlands

Tel number: +31 345 712 070

Email adres: info@goffinmt.com

Website: www.goffinmolculartechnologies.com