



# SMART PARATUBE™

For a cleaner and more efficient way to test for parasites




## Smart Paratube™

MedSchenker™ Smart Paratube™ is a revolutionary solution for parasite egg testing. The Smart Paratube™ is a quick and clean solution for any size lab or testing facility and dramatically improves the testing process by saving 90% of your existing medical waste disposal costs.

The inner tube utilizes a standard mesh filter to collect parasites and extract debris while the outer tube is used for inspection. This allows the lab technician to perform the test without being exposed to any fecal matter.

## Smart Paratube™ (PARA15)

|  |                     |   |
|--|---------------------|---|
|  | <b>Package</b>      | 1000 units per master package   |
|  | <b>Vial</b>         | Screw cap with double tubes   |
|  | <b>Intended use</b> | Parasite egg detection  |
|  | <b>Components</b>   | Outer tube, Inner tube, Outer tube cap, Mesh (1st filter, 2nd filter)               |
|  | <b>Dimensions</b>   | Outer tube: 12.4cm ø 1.8cm<br>Inner tube: 9.5cm ø 1.5cm<br>Mesh (filter): 400-500um |



## ADVANTAGES

- ❄ Dramatically improves the troublesome process of testing for parasites. This revolutionary solution can even be used in a clean environment.
- ❄ Convenient and easy to process. Prepare sample for inspection within 15 minutes.
- ❄ Saves up to 90% of your existing medical waste costs.
- ❄ Cleaner working environment. Perform inspection without being exposed to any fecal matter.
- ❄ Suitable for any size lab or testing facility.

PARA15-0001-Q121



Industrieweg 24C, 4153 BW, Beesd, Nederland



+31 (0)345 712 070



info@goffinmt.com



goffinmoleculartechnologies.com