

URINE COLLECTION SYSTEM

Stackable. Secure. Simple.




Urine Collection System

The MedSchenker® Urine Collection System is a complete system made out of 100% polypropylene and prelabeled.

Our 50ml graduated wide mouth urinalysis cup is stackable, which helps save valuable space and a knurled cap improves gripping while wearing gloves. The 10ml stand-up conical tube is ideal for sediment collection and our CQ-Lock™ screw cap prevents leaks even in pneumatic tubes and centrifuges. For an uncomplicated transfer of the specimen from the urinalysis cup to the tube we provide a 3ml disposable pipette.

Engineered with your needs in mind, the MedSchenker® Urine Collection System is the best choice for your testing needs.

Urine Collection System (UCS-KIT)

	Package	250 units per master package
	Material	Polypropylene
	Cup	50ml Graduated and Stackable
	Tube	10ml with CQ-Lock™ screw cap
	Pipette	3ml Disposable

ADVANTAGES

- 50ml graduated wide mouth urinalysis cup is stackable which saves valuable space.
- 10ml stand-up conical bottom tube is ideal for sediment collection.
- CQ-Lock™ screw cap prevents leaks even in pneumatic tubes and centrifuges.
- Made from 100% polypropylene.
- Urinalysis cup and tube are prelabeled.

UCS-KIT-0001-Q121