

Automated Fully Enclosed qPCR Instrument

Galaxy Neo

Processes **2 samples** simultaneously
Melting Curve Analysis
LIS/HIS Support

1

7 colors Fluorescence channel

2

Built-in Sample Pre-treatment with
ultrasound and Heating Lysis

3

Extraction, Purification, Amplification
and Results Analysis integrated
in one system

4



2 cartridges simultaneously

Ready to use

Pre-loaded reagent Automated
lyophilization technology

1

Fully Enclosed

Contamination free, Biosafety
No PCR lab required

2

Open System

Compatible with third-party PCR kits
for more customized possibilities


3

Compatibility

Extract nucleic acids from complex
biological samples

4

Workflow



Safer · Faster · Easier · More Accurate